

Rhodium Ltd  
4b Sheffield St, Katikati, 3129  
New Zealand  
Email: nzlabs@rhodium.com  
Free Phone: 0800 746 634  
www.rhodium.com

Standard delivery time frame is up to 10 working days\*  
from receiving the impression at our laboratory.  
We will notify you when we have received this.

**1 Dentist/Practice Details**  
Dentist name: \_\_\_\_\_ Date: \_\_\_\_\_  
Practice name: \_\_\_\_\_

**2 Contact and Shipping Details**  
Contact phone number: \_\_\_\_\_  
Shipping address: \_\_\_\_\_  
\_\_\_\_\_


**3 Patient Details**  
Name/Code: \_\_\_\_\_ OVC Placement Date: \_\_\_\_\_

**4 Disinfection confirmation**  
I have disinfected the impression with: \_\_\_\_\_ by: \_\_\_\_\_

**5 Dental Onlay Specifications** (Please circle the tooth number)

U  
R | L  
L

**Profile (Molars only):**  
High Anatomy  Low Anatomy (More conservative)



Cross section of hard layers

**Shade:**

A1 
A2 
A3

**6 Accessories** (1 McDonald Matrix Band, 2 OVC Proximal Wedges and 2 Green Wedges are included complimentary)

Add 2 extra OVC Proximal Wedges:	<input type="checkbox"/> (tick)	\$1
Add 1 extra McDonald Matrix Band:	<input type="checkbox"/> (tick)	\$2
Reusable Hex Sticks:	<input type="text"/> (quantity of packs)	\$29 per pack (contains 2 Hex Sticks)
OVC Contact Instrument:	<input type="text"/> (quantity)	\$39 per instrument
Light Cured Block Out Resin 1.5ml:	<input type="text"/> (quantity)	\$19.95 per tube

Custom made device | Supplied in an unsterilised state

This medical device has been manufactured to satisfy the design characteristics and properties specified by the prescriber for the above referenced patient.  
All prices quoted are GST exclusive.

\*as we use standard shipping services, any delays in shipping are out of our control.

# Impression Requirements



The impression is taken of a un-prepared tooth. It is best taken at the initial consultation. If, at the consultation, the patient is unsure as to which treatment option to take, we recommend taking a “no obligation quick impression” to prevent the patient having to return. This task could be given to a skilled dental assistant to save time.

We recommend dual-arch trays.

Please ensure the following:

- The patient’s teeth must be in centric occlusion.
- Include teeth on either side of the treatment tooth.
- Canines should be in the impression as our lab needs canine-guidance to get the occlusion correct.
- The entire un-prepped crown of the tooth to be treated should be in the impression. Less impression material is required for the opposing arch as only the occlusal surfaces are required.

## Material choice

Any silicone impression material is suitable. We recommend fast set PVS materials such as bite impression materials. Cavex ColorChange Alginate is also acceptable as this is dimensionally stable for 7-9 days. Wrap alginate in damp tissue and place in a Ziplock bag.

If you are concerned about the longevity of alginate, please cast the treatment side with stone-plaster.

Wrap impressions carefully in bubble-wrap to minimise the chance of distortion during shipping to the lab.

---

# Shipping Instructions

Wrap your impression and place it inside appropriate packaging with the completed prescription form.

Shipping address: Rhodium Ltd  
4b Sheffield St, Katikati, 3129  
New Zealand

For terms and conditions visit [www.rhodium.com/tc](http://www.rhodium.com/tc)

**NZ Free Phone: 0800 746 634**